



A novel MR-guided treatment planning tool reduces accumulated critical organ dose in prostate cancer

Ellen H.A. Loeters, Roel G.J. Kierkels
Radiotherapiegroep, Arnhem/Deventer, the Netherlands



Introduction

MR-guided adaptive radiotherapy (MRgART) ensures target coverage and optimal organ at risk (OAR) sparing at each fraction. However, evaluation of accumulated over- or under-dosage over treatment fractions, has not been integrated into the current clinical workflow of the 1.5 T Unity MR-Linac.^[1]

Purpose

- To develop and evaluate a dose accumulation pipeline to estimate and compare the accumulated given dose (D_{Acc}) with the pre-treatment planned dose (D_{PT}), for hypo-fractionated prostate cancer patients;
- Additionally, we aimed to introduce a novel (adaptive) treatment planning workflow for the pre-treatment and adaptive MRgART plans.

Materials/Methods

Clinical IMRT treatment plans of six randomly selected prostate cancer patients, treated with the 1.5 T Unity MR-Linac (Elekta AB, Stockholm), receiving 35 Gy in 5 fractions over 2 weeks, were included. The pre-treatment plans, adaptive plans and contours from the Monaco TPS (TPS_{MON}) (Monaco v5.51.11) were used to evaluate the dose accumulation pipeline in RayStation (TPS_{RAY}) (RayStation v2023-R, RaySearch Laboratories AB, Stockholm). In TPS_{MON} the adaptive planning was based on the adapt-to-shape and adapt-to-position procedure.^[2] Pre-treatment and adaptive plans were also created in TPS_{RAY} and compared with the clinical plans (table 1).

Dose accumulation pipeline

- Hybrid ANACONDA deformable image registrations with controlling ROIs: CTV, bladder and rectum
- TPS_{MON} and TPS_{RAY} fraction doses mapped onto the pre-treatment reference MRI on a voxel-by-voxel level.
- Dosimetric evaluation of D_{Acc} in TPS_{RAY}

RayStation MRgART planning

- Beam modeling for the 1.5 T MR-Linac with crosst and MRI coils
- Deformable CT mapping to T2-weighted MRI for dose calculations
- Final Monte Carlo dose with 1.5 T B-field and 1% MC uncertainty
- Adaptive plans: warm restart of pre-treatment plan (4 x 40 iterations)

Table 1. Intensity modulated MRgART treatment plan optimization settings

| TPS | Min. segment area cm ² | Min. Segment width cm | Min. MU / segment | Max. Number of segments |
|------------|-----------------------------------|-----------------------|-------------------|-------------------------|
| Monaco | 6.0 | 0.7 | 7.0 | 70 |
| RayStation | 4.0 | 0.3 | 4.0 | 70 |

Results

For all patients, the V95% of the CTV was >99.0% in the D_{PT} and the D_{Acc} in both TPSS. A T2-weighted MRI overlaid with a dose difference distribution between the TPS_{MON} and TPS_{RAY} of a typical patient illustrates lower peripheral doses in the TPS_{RAY} plan (Figure 1A). In D_{PT} on average (1 SD), the rectum V28 was 1.5 (1.2) % lower in the TPS_{RAY} plans. The D_{mean} of the evaluated OARs was (range) 0.1 – 6.8 Gy lower in TPS_{RAY} ($p < 0.001$) for all patients compared to TPS_{MON} (Figure 2A). Lower OAR doses may be attributed to a higher degree of modulation within TPS_{RAY} plans. On average, the number of beam segments was 51 (range: 45 – 57) in TPS_{MON} and 69 (68 – 70) in TPS_{RAY} .

On average, the OAR D_{mean} between the D_{PT} and D_{Acc} was similar ($p = 0.26$). However, for one patient, a higher bladder D_{mean} was observed due to inaccurate bladder deformation (Figure 1B and 2B).




Figure 1. (A) Sagittal cross section of a T2-weighted MRI overlaid with a dose difference distribution between the Monaco and RayStation pre-treatment doses (D_{PT}) [red colors indicate lower doses in the RayStation plan], and (B) between the D_{PT} and the accumulated dose. The corresponding dose volume histograms illustrate the target and OAR doses of the (AII) Monaco D_{PT} (dashed lines), the RayStation D_{PT} (solid lines), and the accumulated dose distribution (dotted lines) (BII). Line colors correspond to the figure 2 legend.




Figure 2. Scatter plots showing the relationship between (A) the OARs D_{mean} of the Monaco pre-treatment dose distributions (D_{PT}) and the RayStation D_{PT} and (B) between the OAR D_{mean} of the RayStation D_{PT} and the estimated actually given accumulated dose (D_{Acc}). The arrow points towards the bladder D_{mean} of the patient shown in figure 1b.

Conclusion

- The dose accumulation pipeline provides an estimate of the actually given dose, which was demonstrated for prostate cancer patients treated at the 1.5 T MR-Linac. However, the quality of the deformable registrations should be critically reviewed.
- The novel pre-treatment and adaptive planning class solution in RayStation results in prostate IMRT plans with lower mean OAR doses as compared to the clinical plans created in Monaco.

Acknowledgements

We would like to thank Agnes Angerud (RaySearch Labs), for providing a RayStation version with MR-Linac planning functionalities.

References

- McDonald BA, Zachiu C, Christodouleas J, Naser MA, Ruschin M, Sonke JJ, et al. Dose accumulation for MR-guided adaptive radiotherapy: From practical considerations to state-of-the-art clinical implementation. *Front Oncol* 2023;12(January):1–14.
- Winkel D, Bol GH, Kroon PS, van Asselen B, Hackett SS, Werensteijn-Honingh AM, et al. Adaptive radiotherapy: The Elekta Unity MR-linac concept. *Clin Transl Radiat Oncol* 2019;18:54–9.

Contact

r.kierkels@radiotherapiegroep.nl